


# A SOLUTION TO MEET YOUR AV ACCESS NEEDS

**IN.PACT™ Admiral™**  
Drug-Eluting Balloon



Maintaining access to the arteriovenous circuit is critical to extending the lives of patients. Fortunately advances in drug-eluting balloon therapy provide options to keep the AV access site open.

## **Address AV Access Stenosis with Paclitaxel**

The most common cause of access dysfunction is venous stenosis caused by neointimal hyperplasia. Preclinical studies show paclitaxel is beneficial when treating AV stenosis.

## **Extend Time to Reintervention**

Physician initiated studies show promising clinical outcomes with IN.PACT Admiral relative to other AV access therapies.

## **Treat with Ease - New 40 cm Catheter**

Short 40 & 80 cm catheters allow for shorter wires and simplified handling.

# IN.PACT™ Admiral™ Drug-Eluting Balloon

Ref. Number 40 cm	Ref. Number 80 cm	Ref. Number 130 cm	Balloon (mm)	Balloon Length (mm)	Recom. Introducer sheath (F)	RBP
SBI 040 040 04P	SBI 040 040 08P	SBI 040 040 13P	4	40	5	14
SBI 040 060 04P	SBI 040 060 08P	SBI 040 060 13P	4	60	5	14
SBI 040 080 04P	SBI 040 080 08P	SBI 040 080 13P	4	80	5	14
SBI 050 040 04P	SBI 050 040 08P	SBI 050 040 13P	5	40	6	14
SBI 050 060 04P	SBI 050 060 08P	SBI 050 060 13P	5	60	6	14
SBI 050 080 04P	SBI 050 080 08P	SBI 050 080 13P	5	80	6	14
SBI 060 040 04P	SBI 060 040 08P	SBI 060 040 13P	6	40	6	14
SBI 060 060 04P	SBI 060 060 08P	SBI 060 060 13P	6	60	6	14
SBI 060 080 04P	SBI 060 080 08P	SBI 060 080 13P	6	80	6	14
SBI 070 040 04P	SBI 070 040 08P	SBI 070 040 13P	7	40	7	14
SBI 070 060 04P	SBI 070 060 08P	SBI 070 060 13P	7	60	7	14
SBI 070 080 04P	SBI 070 080 08P	SBI 070 080 13P	7	80	7	14

## Medtronic

**Invatec S.p.A.**  
Via Martiri della Libertà 7  
25030 Roncadelle (BS)  
Italy  
Tel: +39.030.2589311

**Medtronic**  
710 Medtronic Parkway NE  
Minneapolis, MN 55432  
USA  
Tel: +1.763.514.4000

**Medtronic International  
Trading Sàrl**  
Route du Molliau 31  
CH-1131 Tolochenaz  
Switzerland  
Tel: +41.21.802.7000

**Medtronic of Canada Ltd.**  
99 Hereford Street  
Brampton, Ontario L6Y 0R3  
Canada  
Tel: +1.905.460.3800

**Medtronic Latin America**  
3750 NW 87th Avenue, Suite 700  
Miami, FL 33178  
USA  
Tel: +1.786.709.4200

**Medtronic International Ltd.**  
49 Changi South Avenue 2  
Singapore 486056  
Tel: +65.6436.5000

**Medtronic Australasia Pty Ltd.**  
97 Waterloo Road  
North Ryde, NSW 2113  
Australia  
Tel: +61.29857.9000

**Medtronic New Zealand**  
Unit N16, Mezzanine Level 5  
Gloucester Park Road  
Onehunga  
Auckland  
New Zealand

**Medtronic Korea Co., Ltd.**  
5F, Sajjo Building  
1001 Daechi-dong, Kangnam-ku  
Seoul, 135-280  
Korea  
Tel: +82.2.3404.3600

**Medtronic META FZ-LLC**  
Office Park, Block D, 2nd Floor  
PO Box 500638  
Dubai Internet City | Dubai,  
United Arab Emirates  
Tel: +971.4.818.2666

[medtronic.com/apv](http://medtronic.com/apv)